

Pricelist of Oligonucleotides on Microplate

Prices in **EUR** (VAT excluded)
pricelist updated:
4.1.2010

Basic Pricelist

item	price	
Set of 96 oligonucleotides (up to 30 bases, incl.) on 96-well-microplate	680	price per microplate
Set of 96 oligonucleotides (up to 40 bases, incl.) on 96-well-microplate	870	price per microplate
Set of 96 oligonucleotides (up to 50 bases, incl.) on 96-well-microplate	1 060	price per microplate
Set of 96 oligonucleotides (up to 60 bases, incl.) on 96-well-microplate	1 250	price per microplate

surcharges to the price of an microplate

item	price	
Surcharge for unification of molar amounts of all oligos (equimolar microplate)	76	surcharge per microplate
Surcharge for 5'-amino modification	230	surcharge per microplate
Surcharge for inosin (for every 100 inosins)	110	surcharge for every 100 inosins

The output is a standard 96-well microplate with dry oligonucleotides whose sequences and locations on the array have been defined by customers when ordered. The volume of the well (2 mL) enables the oligo to be ordinarily dissolved.

The molar amounts of the oligos in the wells are either unadjusted (they differ among the wells), or unified to one equimolar value (then it is possible do dissolve all the oligos by even volume of water or buffer).

