

# GEMINI™ probes

GEMINI™ probes are a type of real-time PCR probe developed in GENERI BIOTECH. The technology of GEMINI is based on placing internal quenchers inside the probe sequence that results in a stronger quenching and a lower background signal. The probes are HPLC purified (no additional fee).

**Amounts available:** 5 nmol and 10 nmol

Product	Amount of probe	
	5 nmol	10 nmol
GEMINI™ FAM	1507-005	1507-010
GEMINI™ HEX	1505-005	1505-010

**Attention:** probes containing a fluorophore are recommended to be stored in the dark.

*In the table are stated catalogue numbers of modifications*